

MICROCOPY RESOLUTION TEST CHART
NATIONAL BUREAU OF STANDARDS-1963-A

OTIC FILE COPY

REPORT NUMBER 4

MEASURING SLEEP BY WRIST ACTIGRAPH

FINAL REPORT

Daniel F. Kripke, John B. Webster, Daniel J. Mullaney, Sam Messin, and Paul Fleck

June 1, 1981

Supported By

U.S. ARMY MEDICAL RESEARCH AND DEVELOPMENT COMMAND

Fort Detrick, Frederick, Maryland 21701

Contract No. DAMD-17-78-C-8040
University of California
La Jolla, California 92093



Approved for public release; distribution unlimited

The findings in this report are not to be construed as an official Department of the Army position unless so designated by other authorized documents.

83 02 07 011

TABLE OF CONTENTS

	page
ABSTRACT	2
INTRODUCTION	2
EXPERIMENT 1: SCORING SLEEP FROM WRIST ACTIVITY RECORDINGS	3
EXPERIMENT 2: TRANSDUCER OPTIMIZATION	4
FIGURES 1 AND 2	5-7
FIGURE 3	8
EXPERIMENT 3: TRANSDUCER ORIENTATION	9
EXPERIMENT 4: TRANSDUCER PLACEMENT	9
FIGURE 4	10
FIGURE 5	11
EXPERIMENT 5: DIGITAL PREPROCESSING	12
FIGURE 6	13
EXPERIMENT 6: SLEEP RECOGNITION	14
FIGURE 7	15
TABLE 1	16
TABLE 2	19
A WEARABLE SLEEP MONITOR SYSTEM	20
FIGURE 8	21
FIGURE 9	23
SLEEP RECOGNITION SOFTWARE	24
ASSESSMENT OF VALIDATION	25
TABLE 3	26
FIGURE 10	27
TABLE 4	28
HARDWARE SPECIFICATIONS OF A WRIST-WEARABLE SLEEP MONITOR	3 0
SOFTWARE SPECIFICATION	31
SUMMARY	31
ACKNOWLEDGMENT	33
REFERENCES	34
ACKNOWLEDGMENT REFERENCES APPENDIX 1 Accompany of the second seco	35-44
APPENDIX 1	
y. Je le le	
July or Codes	
Distribution Codes Notation Special	
Dieta Ispirit sun)
1 0150	~
· / / / /	

ABSTRACT

Future warfare may be exceptionally intense and brief, offering personnel little time for sleep. Fatigue could be a crucial factor in combat performance, especially if troops are airlifted across time zones and subjected to jet lag. Therefore, monitoring sleep loss and facilitating adequate sleep are crucial aspects of military medical planning. To monitor sleep in field conditions, a new technology is needed.

People move constantly when awake but little when asleep. Electronic recognition of activity can thus be used to monitor and determine sleep/wake states. This Final Report summarizes the research and development of a wearable sleep monitoring system. Wrist activity is measured with a piezo-ceramic transducer, monitored and stored by a microprocessor, then transferred to a larger computer for automatic sleep scoring. In prospective validation trials, the automatic measurement of sleep by a prototype device correlated r=0.97 with BEG sleep scoring. Full technical specifications are presented for construction of field-deployable sleep monitors which could be worn entirely on the wrist. Deployment of such devices would permit operational objective measurement of sleep loss among our troops.

INTRODUCTION

Sleep loss and combat fatigue are increasing concerns to the modern army. Lapses of attention may be more serious for the soldier who operates smart weapons, vehicles, radar, etc, than they were for the traditional infantryman who carried a pack and dug foxholes. Any future war is likely to be extremely brief and intense, and soldiers using technologically sophisticated weaponry will have little time for sleep. American troops may have to enter combat immediately after being airlifted to remote parts of the world, and plans must be developed to minimize the effects of jet lag on personnel performance.

In order to assess the effects of sleep loss, sleep schedule changes, jet lag, etc, upon operational performance of troops in actual field conditions, a practical method of quantifying sleep is needed which is cost-effective and reliable. Traditional methods of measuring sleep through EEG-EOG-EMG recordings are impractical for actual or simulated combat settings, and subjective monitoring is unreliable (Carskadon et al, 1976). Electrophysiological and observational methods for measuring sleep are also costly, and considerable time is necessary to score EEG records.

In the three years of our contract (DAMD 17-78-C-8040), we have investigated the feasibility of recognizing sleep from measures of wrist activity. We have developed instruments and procedures which allow us to quantify sleep with accuracy comparable to EEG scoring, but with much less expense and inconvenience. In this Final Report, we review and update the research and development reported in prevous Annual Reports and include specifications for building sleep monitors suitable for field deployment. The

next step would be to actually build and deploy miniaturized monitors in operational settings. The technical specifications necessary to construct such monitors have been developed, and the monitors could be built today with a suitable engineering investment.

EXPERIMENT 1: SCORING SLEEP FROM WRIST ACTIVITY RECORDINGS

Kupfer et al (1972) and Foster et al (1972a, 1972b) described the use of a telemetric activity recording system for quantifying sleep. Kripke et al (1977) developed a more flexible system for recording wrist activity using a piezo-ceramic activity transducer and Medilog* cassette recorder. Initial pilot results established the value of this analog recording device for quantifying sleep. Kripke et al (1978), using five healthy subjects, obtained a correlation of 0.98 between sleep duration determined from wrist activity and from the EEG. The first experiment to be described here demonstrated the reliability of analog wrist activity recordings for quantifying sleep with a much larger number of recordings and a more diverse group of subjects, including hospital patients with varying degrees and types of insomnia.

Method

The activity transducer was constructed by soldering a small steel nut onto a 5 mm length of spring-like polygraph pen cleaning wire, and clamping the other end of the wire against a piezo-ceramic element (Kripke et al, 1978). The transducer was packaged in a small acrylic box mounted on a watchband. Transducer output was recorded on one channel of a four-channel Medilog cassette recorder. EEG, EOG, and EMG were simultaneously recorded on the other three channels, with timing marks superimposed on the EMG channel. In all, 102 overnight recordings were obtained from 85 subjects, including 39 recordings from 32 hospital patients.

The EEG-EOG-EMG signals (with the time code) were replayed from the cassette to a polygraph at 15 mm/second. Independently, the amslog activity recording for the same segment was replayed (with the time code) onto a separate chart at 32 mm/minute. The EEG-EOG-EMG records were scored by one investigator, while activity records were scored independently by a second investigator, blind to the EEG results.

Results and Discussion

Figure 1 shows a representative analog activity recording. Three different sections of the same record (separated by dots) illustrate sleep onset, midsleep, and sleep offset, respectively. Straight lines indicate inactivity while modulations of the line indicate activity. The scorer's sleep/wake determination ('S' or 'W') is noted about the tracings. Periods of sleep and wake can be readily identified by the relative amounts of activity recorded.

The number of minutes for which the activity scorer and the EEG scorer agreed in assigning a score of "wake" or "sleep" was determined for each record. For the entire sample of 102 records, the two scorers agreed on 94.5% of all minutes. Agreement for the 63 non-patient records was 96.3%. The correlation coefficient between the two estimates of Total Sleep Time (TST) was r=0.90,

^{*}Ambulatory Monitoring, Inc., 731 Saw Mill River Rd., Ardsley, N.Y. 10502.

indicating that most of the variance in inter-subject TST scored by one method was represented by the other method. The activity scorer's estimate of TST averaged 15.33 minutes greater per night (t=3.82, p<0.0001). This trend was related to a tendency to underestimate sleep onset latency in activity records. The activity rater scored sleep onset first in 54 records, while the EEG scorer scored sleep onset first in 34 records. The two scorers agreed exactly on sleep onset in the remaining 14 records. On average, sleep was scored four minutes earlier in activity than EEG records.

Figure 2 compares activity scoring (top line) with EEG scoring (bottom line) for ten records (including seven patient records). Thin lines indicate wake. Disagreement is indicated by a dark area below each pair of lines. Typical were the relatively good agreements between sleep/wake scoring for minutes scored from EEG and from activity recordings, as well as the high correlation between TST estimated by both methods and the small absolute error in TST per night. This confirmed the earlier findings and encouraged further development of a sleep monitor system based on activity recording (Mullaney et al, 1980).

EXPERIMENT 2: TRANSDUCER OPTIMIZATION

Although the piezo-ceramic activity transducer developed by Kripke et al (1977) produced excellent recordings and allowed highly accurate sleep/wake scoring in Experiment 1, the possibility that other available activity transducers might be superior was considered. Experiment 2 compared the piezo-ceramic transducer to a commercially available tilt-switch motion transducer* and to a sensitive accelerometer.**

Method

The tilt-switch transducer and a piezo-ceramic transducer were mounted in a 3.7 x 3.5 x 5.6 cm acrylic box and connected to two channels of a Medilog recorder. A 1.35 V battery and a resistive voltage divider were necessary to match the output of the tilt-switches to the input requirements of the recorder. Four subjects were the pair of transducers on their wrist for a total of six nights. The two channels were played back simultaneously onto polygraph paper at 32 mm/minute. In further experiments, a second piezo-ceramic transducer in a similar box was firmly taped to the accelerometer. The two devices were connected to two channels of the polygraph and were worn by two subjects for several hours each. The duration of these recordings was limited since the accelerometer had to be directly connected to the polygraph.

Results and Discussion


Figure 3 shows a representative example of the performance of the piezo-ceramic and tilt-switch activity transducers. The piezo-ceramic transducer detected movements repeatedly at times when the tilt-switch transducer did not. Apparently, many movements did not change the orientation of the tilt switch array sufficiently to throw any of the switches. In our entire sample, there were no examples where the tilt-switch transducer detected activity not recorded by the piezo-ceramic transducer.

^{*}Vitalog Corporation, 1058 California Avenue, Palo Alto, CA 94306. **Grass Instruments, Model SPA 1.

Figure 1. Representative analog activity record visually scored in Experiment 1. Activity is seen as modulations of the line, while straight lines indicate inactivity. The scorer's marks ('S' and 'W') are seen on the record. The first three lines show sleep onset, the middle line mid-sleep, and the last three lines show awakening.

Figure 2. Comparison of the scoring of 10 records from analog activity data (top line) and EEG-EOG-EMG data (bottom line). Thick lines indicate sleep, thin lines wake. Disagreement is indicated by a dark space beneath each pair of lines. Seven of the ten records are from patients.

AND THE REPORT OF THE PROPERTY OF THE PROPERTY OF



EEG-EOG-EMG Activity 7 63




Figure 3. Representative polygraph record of activity recorded simultaneously from piezo-ceramic activity transducer (Channel 1, top) and Vitalog tilt-switch activity transducer (Channel 2). Although the piezo-ceramic transducer measures activity throughout the record, the tilt-switch fails to detect much of this activity.

A representative accelerometer recording is presented in Figure 4. The piezo-ceramic transducer detected all movements detected by the precision accelerometer, but the accelerometer was not sensitive to all movements, presumably due to its uniaxial design.

Thus, the piezo-ceramic transducer compared favorably to alternative designs, yet it is extremely compact, inexpensive and easy to manufacture.

EXPERIMENT 3: TRANSDUCER ORIENTATION

Although the weight is mounted off-center on the spring wire, causing the piezo-ceramic transducer to be excited by accelerations or rotations in any axis, the output voltage was not equalized for each axis. In Experiment 3, transducers were mounted in different axes on the wrist to determine optimal transducer orientation.

Method

Three piezo-ceramic transducers were mounted at three perpendicular orientations within a single 3.7 \times 3.5 \times 5.6 cm box connected to three channels of a Medilog recorder. Six subjects were this device on their right wrists for a total of ten nights. The three channels were replayed simultaneously onto polygraph paper at a rate of 32 mm/minute.

Results and Discussion

An example of a recording from the three transducers is presented in Figure 5. This example is representative of the entire sample. Although the recorded activity analog from one transducer was frequently somewhat larger than that from another, there were virtually no instances where movements were not registered by all three transducers, and no clear superiority of orientation could be determined. It was concluded that orientation is not critical for this transducer.

EXPERIMENT 4: TRANSDUCER PLACEMENT

The decision to mount activity transducers on the dominant wrist in previous work was based on the a priori assumption that this location would allow the most sensitive detection of activity. To determine whether mounting on the non-dominant wrist or another area might detect more movements, we simultaneously recorded from both wrists, an ankle, and the forehead. The forehead was chosen with the supposition that all truncal hovements would be accompanied by head movement, but respiratory artifacts might be less pronounced on the head than on any site on the body trunk.

Method

Four transdur s were mounted in separate boxes (3.5 x 4.4 x 1.7 cm) and connected withe our channels of a Medilog recorder. Mine subjects completed




Figure 4. Representative polygraph record of piezo-ceramic transducer (Channel 1, top) and precision accelerometer (Channel 2). Piezo-ceramic transducer detects all activity measured by accelerometer, and some activity missed by accelerometer.




Figure 5. Representative polygraph record of 3 piezo-ceramic transducers mounted at the 3 major axes within a single box and worn on the wrist. Although differences in amplitude occur between the 3 channels, there are no failures to detect activity in any channel.

22 overnight recordings with the transducers worn simultaneously on each wrist, the forehead and an ankle. The assignment of transducers to locations was counterbalanced to control for possible variations in the sensitivity of individual transducers, and the order of polygraphic presentation was also counterbalanced to control for rater bias. The four channels were replayed simultaneously to four channels of the polygraph at 32 mm/minute. The records were then examined by two experienced raters (blind to transducer location) to score each transducer recording for sensitivity to movement and utility for sleep/wake scoring.

Results and Discussion

Figure 6 is part of a typical four-channel recording, showing more movements in the wrists than in either the head or ankle. The raters felt the contrasts between the two wrists were usually very small and of little importance in sleep scoring, and either wrist was preferable to the head or ankle. There was evidence that the left (non-dominant) wrist produced somewhat more activity. The practice of recording wrist activity is supported by these results. The choice of wrists seems unimportant, but the more comfortable non-dominant wrist is probably the better choice.

EXPERIMENT 5: DIGITAL PREPROCESSING

In Experiments 1 to 4, analog activity recordings were collected on cassette tape, replayed to a polygraph, and then scored visually. With this analog system, huge amounts of largely redundant data were collected, much time was expended, both in replaying and scoring recordings, and delicate instruments and trained technicians were required. A much faster and more efficient system would store a measure of activity in a digital memory and make it accessible to a computer for automatic scoring. The first problem in implementing such a system was to find an effective method of expressing a period of continuous activity recording as a single digital value. This experiment compared ten approaches for data compression and digital storage.

Method

Recorded analog activity data were replayed on one channel of a polygraph, and, at the same time, sent through an analog-to-digital converter (ADC) to a Hewlett-Packard 2100 computer. The computer was programmed with ten different data compression algorithms which calculated ten different digital values for each 2-second data epoch. The computer also generated a time code each minute which was written to a second polygraph channel.

The ADC converted the analog signal to a digital value 240 times a second. The 240 Hz conversion rate was chosen to be exactly four times the frequency of 60 Hz electrical noise occasionally recorded (such as that produced by an electric blanket or clock near the bed during the night.) It will be shown that the sum of every four conversions at 240 Hz cancels 60 Hz interference.

A file of seven overnight wrist activity recordings was digitized and transformed into ten alternative digital files corresponding to the ten




Figure 6. Representative polygraph record of activity detected by piezo-ceramic transducers mounted on the right wrist (Channel 1), left wrist (Channel 2), right ankle (Channel 3), and forehead (Channel 4). Bilateral wrist activity frequently occurs in the absence of head or ankle motion.

algorithms. Portions of the digitized records were displayed on a plotter to compare the ten transformations during different forms of activity. Each of the ten alternative digital files was then scored automatically by the computer using a simple rule in which a minute was scored 'wake' if x of the 30 two-second epoch scores exceeded a threshold of y. Sleep/wake scores for a range of x and y values were compared to the scores assigned by an experienced scorer to the polygraph activity records.

Results and Discussion

Figure 7 shows the plotter display and polygraph write-out of a five-minute portion of a record contaminated with 60 Hz noise from an electric blanket. The ten horizontal traces on the plot represent the ten digital transformations of the analog activity displayed on the polygraph. The vertical lines on the plot separate minutes, which are also marked and labelled with a binary code on the polygraph paper. Of particular interest in this figure is the contrast between traces 1, 3, 5, and 7 and traces 2, 4, 6, and 8 during periods of electric blanket noise. Even-numbered transformations squared conversions prior to summation of each 4 measurements, so all sums were positive and cancellation of 60 Hz noise did not occur. Trnasformations 1, 3, 5, and 7 summed conversions prior to squaring so cancellation of noise can and did occur. The absence of noise in these latter traces indicates the effectiveness of the simple digital filtering technique.

The potential of each of the ten methods of digital preprocessing was evaluated by calculating the maximum agreement between sleep/wake status computed automatically from each digital file and scored from the polygraph activity record. Table 1 presents the rank order of maximum agreement for each record. Despite some variability, the indication from these data is that transformation 5 was superior to the others. (Transformations 1, 2, and 10 were not tested completely since they were judged inadequate after preliminary evaluations). Transformation 5 summed the absolute difference between the activity score for each 1/60th second and the mean of the five preceding and five subsequent scores.

To summarize the preprocessing algorithm found optimal in this experiment, the analog output of the activity transducer was digitized at 240 Hz, and every four digital activity values was summed to produce a value free of any 60 Hz component. Each of the resulting 120 scores (y) in every two-second epoch were then converted to 120 difference scores, y(diff)=10 * y(1) - [y(1-5) + y(1-4) + y(1-3) + y(1-2) + y(1-1) + y(1+1) + y(1+2) + y(1+3) + y(1+4) + y(1+5)], where y(1) is the current 60 Hz sum, <math>y(i-1) the preceding sum, y(i+1) the subsequent sum, etc. The sum of absolute values of these 120 difference scores was the final quantitative expression of activity over a two-second data epoch. Although these algorithms were tested by a retrospective procedure not strictly comparable to prospective scoring, it was encouraging that the median percent agreement obtained for the best algorithm was 91%.

EXPERIMENT 6: SLEEP RECOGNITION

The sleep recognition procedure used in Experiment 5 was an extremely




Figure 7. Five-minute portion of analog-polygraph record (top) showing time code and activity mixed with 60 Hz interference from an electric blanket, and plotter display (bottom) showing the ten digital transformations of the same 5-minute segment of activity data. Transformations 1, 3, 5 and 7 effectively filter the 60 Hz noise. A and B are 60 Hz interference.

Table 1.

Rank ordering of agreement scores for 7 of the 10 digital transformations of activity data for each of the 7 records analyzed in Experiment 5.

		Transform Number						
		<u>3</u>	4	<u>5</u>	<u>6</u>	<u>7</u>	<u>8</u>	<u>9</u>
Record Number	1	5	2	6	3	7	1	4
	2	7	2	1	2	2	2	2
	3	6	5	1	3	2	4	7
	4	7	4	1	5	5	3	1
	5	2	2	6	5	6	2	1
	6	2	2	1	2	6	2	7
	7	3	4	1	4	1	4	7
Mean		4.6	2.9	2.4	3.4	6.6	2.6	4.1
Median		5	2	1	3	5	2	4

simple procedure designed to compare preprocessing strategies. Having established an effective preprocessing approach, the next priority was to refine the sleep recognition algorithm for more reliable scoring of sleep.

Method

EEG, EOG, and EMG were recorded from subjects during both wake and sleep, and written onto polygraph paper with a time code. These polygraphic records were then scored for sleep and wake. A wrist activity transducer signal was sampled by the ADC at a conversion rate of 240 Hz. The analog data were preprocessed and stored as described in Experiment 5, but only the optimal transformation (number 5) was used for data analyses. A total of 20 records (13,488 minutes) was analyzed.

Development of the sleep recognition algorithm began with expressions incorporating a weighted sum of combinations of the digital data with potential for discriminating sleep from wake. A minute was scored 'wake' if the weighted sum exceeded a threshold value. For each combination of weights, the agreement between the automatic scoring and EEG scoring was computed, and used as a retrospective measure of the effectiveness of the weighting. The computer varied the weighting of one term at a time iteratively, searching for the combination of weights which produced the highest agreement. Many hundreds of hours of minicomputer time were used to locate the optimal weighting functions.

Seventeen of the 20 records were used in this algorithm development phase. The remaining three records were scored prospectively, ie, each of the three records was scored individually with the single weighting and scale factor found optimal in the development phase. In this prospective test, the laboratory computer simulated the actual deployment of a miniaturized automatic sleep scoring system.

Results and Discussion

One result of the analysis of potential discriminators of sleep and wake was that the most discriminating measure of activity in a minute proved to be the value of the greatest two-second epoch activity amplitude in that minute. For sleep/wake scoring, once the amplitude of the maximal two-second epoch had been entered into the algorithm, no data representing other 2-second epochs from the minute contributed to the discrimination. This unexpected result is extremely fortunate since it allows an entire minute of activity to be summarized by a single value, saving memory space and analysis time.

The optimal sleep recognition algorithm reached after analysis of the 17 records was:

D = 0.025 * [.15y(1-4) + .15y(1-3) + .15y(1-2) + .08y(1-1) + .12y(1) + .12y(1+1) + .13y(1+2)]

where y represented the value of the maximal two-second epoch score in the current minute, y(1-1) that maximal value for the previous minute, y(1+1) for the subsequent minute, etc. If D>1.0, the minute was scored 'wake', otherwise 'sleep'. The agreement between sleep/wake scored automatically with this algorithm and scored from EEG records was 95.56%. Agreement scores and the

proportion of the record scored as sleep by EEG and by the automatic algorithm for each individual record are shown in Table 2.

The ability of this algorithm to score sleep/wake prospectively was tested with the remaining three records. For these records, only the single expression found optimal in the algorithm development phase was chosen prospectively to automatically score sleep/wake. Agreement with EEG scoring and the proportion of the record scored as sleep by EEG and activity are also shown in Table 2. Overall agreement of these three records with EEG scoring was 96.02%.

Since the records analyzed in this experiment varied greatly in length, and included daytime naps as well as overnight sleep, sleep duration was expressed as the proportion of the record scored as sleep (% Sleep). Table 2 includes % Sleep for each record as scored by EEG and as scored automatically from activity data. The automatic scoring method overestimated % Sleep by an average of 1.89% and 2.33% for the retrospective and prospective sets of records. This corresponds to an overestimation of 12.17 and 19.33 minutes per record, or 27.09 and 33.60 minutes per 24-hour day. Correlation coefficients between % Sleep as scored by activity and EEG were r=0.9889 retrospectively and r=0.9982 prospectively. Automatic activity estimates of sleep percentage thus conformed closely with EEG estimates.

The primary automatic scoring algorithm developed here as well as the visual scoring of analog activity records in Experiment 1 tended to overestimate sleep duration as scored by EEG. This was due in part to the fact that sleep onset is delayed somewhat after all activity ceases. Other phenomena causing overestimation of sleep time from activity were recognized from inspection of error listings. For example, a few minutes of inactivity surrounded by many minutes of activity usually did not indicate sleep, although it was sometimes scored as sleep by the primary scoring algorithm. These considerations led to the development of a final stage of the sleep recognition algorithm in which corrections of this nature were made.

In this final stage, the sleep/wake score for each minute as determined by the automatic scoring algorithm was subjected to additional corrections. These corrections took the form "after at least x minutes scored wake, the first y minutes scored sleep are rescored wake" and "y or less minutes scored sleep surrounded by at least x minutes scored wake (before and after) are rescored wake." Substituting a range of x and y values into these expressions and rescoring the records scored retrospectively in Experiment 6, a combination of rules emerged which increased still further the similarity between EEG and activity scoring. These rules were: 1) after at least 4 minutes scored wake, the first minute scored sleep by the primary algorithm is rescored wake; 2) after at least 10 minutes scored wake, the first 3 minutes scored sleep by the primary algorithm are re-scored wake; 3) after at least 15 minutes scored wake by the primary algorithm, the first 4 minutes scored sleep are rescored wake: 4) 6 minutes or less scored sleep by the primary algorithm, surrounded by at least 10 minutes before and after scored wake are rescored wake; and 5) 10 minutes or less scored sleep by the primary algorithm, surrounded by at least 20 minutes before and after scored wake, are rescored wake.

This rescoring procedure increased the agreement of the 17 records scored

Table 2.

Record duration, proportion of the record for which EEG and automatically scored activity scores agree, proportion of the record scored as sleep by each technique, and the correlation coefficient (r) between the proportion scored sleep by the two methods. Activity data was digitized by the laboratory computer. Retrospective records (1-17) were those used to derive the scoring parameters, prospective records (18-20) used the parameters derived for the previous records.

Recor	Recording Dur- ation (mins)	Z Agreement	% Sleep (EEG)	% Sleep (Act.)
1	353	97.17	0.00	2.83
2	574	96.17	45.47	46.86
3	632	95.89	56.80	59.65
4	903	88.82	29.57	34.33
5	660	97.42	54.55	56.21
6	798	95.36	41.73	44.11
7	845	96.80	37.99	40.95
8	1129	9 6.19	30.65	31.62
9	644	96.89	50.47	51.40
10	553	88.79	56.06	48.10
11	371	96.77	92.72	95.96
12	52 7	96.96	14.99	16.13
13	226	89.82	93.81	93.36
14	673	91.98	50.67	45.62
15	593	95.78	56.49	59.36
16	829	91.19	40.17	47.77
17	692	94.08	<u>15.17</u>	20.18
	Total			
1	Retro- 11902 spec-	94.46	42.09	43.98
•	tive		r=0.988	19
18	369	93.50	70.19	76.69
19	846	93.62	28.84	31.68
20	<u>1271</u>	<u>98.35</u>	36.82	37.69
	Total Pros- 2486	96.02	39.06	41.39
3	pecti ve		r=0.998	32

retrospectively from 94.46 to 94.74. The difference in % Sleep was reduced from 1.89 to 0.81, corresponding to 5.24 minutes per record, or 11.65 minutes per 24-hour day. The correlation coefficient decreased negligibly from .9889 to .9868.

A further test conducted with these data sought to determine the resolution in the stored data necessary to achieve these levels of accuracy. The digital activity value was originally stored on disc as a 16-bit value, ie, a number in the range of 0 to 32767. To investigate the bit resolution requirement, the sleep recognition program was repeated with the same data, but the bit resolution was reduced by dividing the data by powers of 2. There was no decrease in agreement with 4-bit data (0-15), and a decrease on only 0.1% was found with 3-bit data (0-7). This surprizing result is important, since it means that data for more minutes can be stored in a given memory capacity. With proper optimization of scaling, 4-bits per minute may be sufficient for reliable sleep/wake scoring. Unfortunately, we have not achieved sufficient experience with the appropriate scale factors to specify which four bits per minute would be most suitable in practice. Certainly 8 bits per minute would be sufficient.

A WEARABLE SLEEP MONITOR SYSTEM

The results summarized to this point were detailed in our previous Annual Reports and have established the feasibility of a sleep monitor system incorporating an activity transducer, digital preprocessor and recorder, and sleep-recognition software. Such a system is capable of quantifying sleep with reasonable accuracy while avoiding the impracticality and expense of polygraphic recording and scoring. Here we describe the implementation of this sleep monitor system in a wearable form.


Digital preprocessing converts the continuous analog signal from the activity transducer to a single digital value for each minute and stores the value in digital memory. To achieve a wearable system, these steps must be accomplished by a digital recorder with the capabilities of a computer, small enough and light enough to be worn comfortably. The monitor must function for days on a battery pack which must also be small and light. To validate these concepts, a wearable monitor was built and tested. Although the prototype was not sufficiently miniaturized for field deployment, it was indeed worn by our test subjects.

To allow validation testing, the Vitalog Corporation* built a package to our specifications incorporating an IM6100 microprocessor, IM6101 programmable interface element (PIE), NSC 809 8-channel analog-to-digital converter (ADC), 6K x 12 random access memory (RAM), crystal clock and a LED indicator light. The entire prototype package, including batteries, is enclosed in a plastic case 15 x 9 x 5-1/2 cm and weighs 480 g. It is about the size of a 35 mm camera and was worn by subjects on a belt around the waist. We built an external activity transducer and connected it to the ADC through a cable. A plug was also provided through which the memory could be accessed by an external Apple II computer. All electronic components were CMOS for minimal power consumption. A block diagram is included (Fig 8).

^{*}Vitalog Corporation, 1058 California Avenue, Palo Alto, CA 94306.

Figure 8.

version, and store the results in RAM each minute. An external Apple microcomputer can be the ADC and RAM is via data bus and control lines organized by the PIE, which also handles Activity (and illumination) connected to the microprocessor to read the data from RAM for permanent storage and sleep generate interrupts which cause the program to activate the ADC, process the digital conscoring. Time information is transmitted through a LED. Communication between the MPU, is detected by the external transducer and transmitted to the ADC. 240 Hz clock pulses Communication with the Apple 1s Block diagram of wearable components of sleep monitor system. interrupt signals, the wait state, and LED output. handled by an EPROM monitor.



In the following sections, the major components of our prototype and their functions are described in more detail. The actual integrated circuit types and manufacturers are listed only for illustration, as similar systems could be designed with CMOS microprocessors of other manufacture.

IM6100 Microprocessor

Overall control of data collection, processing, and storage was accomplished by the IM6100 microprocessor. The IM6100 is a fully programmable 12-bit microprocessor which is functionally equivalent to a PDP-8/e minicomputer, and it recognizes the PDP-8e instruction set. Internal state timing is generated by a crystal which also produces a 240 Hz interrupt signal. A simple accumulator oriented instruction set proved adequate. An 8-bit microprocessor would have been adequate.

An extremely useful feature of the IM6100 is a wait state, in which program execution is suspended but not halted while power consumption is greatly reduced. Battery drain for the entire package (exclusive of the external transducer which contains its own battery) was measured at 8.5 mA (6V) with the processor running continuously, compared with 3.4 mA when halted. By entering the wait state after processing each ADC conversion, battery drain could be reduced to 5.2 mA. At a drain of 5.2 mA, the battery charge was sufficient to operate the monitor program continuously for 180 hours. By entering the halt mode after shorter recording session, the data could be preserved in the digital memory for even longer periods.

IM6101 Programmable Interface Element

The IM6101 PIE handles all communication between the IM6100, the ADC, and LED indicator light, and an external computer, and controls the wait state. Designed for interfacing peripheral devices to the IM6100, the IM6101 allows efficient interaction with these different classes of I/O without the neel for additional external logic, and greatly facilitates interrupt processing. In our reslization, the LED indicator was used to signal time codes by which segments of EEG records could be synchronized with the digital scoring results.


Analog to Digital Converter (ADC)

The ADC converts a voltage in the range of 0 to +5V to an 8-bit digital value (0-255). Up to eight channels of analog data could be sampled by the ADC multiplexer under program control.

External Activity Transducer and Amplifier

As noted above, the ADC converts voltage in the range of 0 to +5V with a resolution of 8 bits. Thus, ground level is converted to '.0', and +5V to '.255'. Maximum resolution is obtained when input voltage covers as much of this range as possible. Since the piezo-ceramic transducer generates only a few milivolts, and alternately produces positive and negative voltage, amplification and level-shifting circuitry powered by a small 5.6V battery was necessary to match the requirements of the ADC. The circuit shown in Figure 9 amplifies the signal from the piezo-ceramic element and shifts it positively so that it alternates

SCHEMATIC DIAGRAM OF EXTERNAL ACTIVITY TRANSDUCER, PHOTOCELL AND LEVEL-MATCHING AMPLIFICATION CIRCUITRY FIGURE 9.



about the level of a center tap in the battery (nominally 2.8V). The signal is amplified so that it covers the full range of 0 to +5V when the transducer is shaken vigorously. This amplified signal is led to the ADC along with the center tap voltage, which serves as a reference. The difference between the signal and reference voltage represents the signed amplitude of the transducer output. A photocell is also mounted in the external transducer package and connected through a resistive bridge to a third channel of the ADC. The voltage level of this photocell circuit serves as a measure of illumination. This can be used as an objective measure of bedtime (defined by "lights out") and can provide some information on the subject's environment.

The Monitor Program

The monitor program (Appendix 1) is loaded from an Apple II microcomputer. After clearing the data memory, the program causes the processor to enter the power-down wait state. On each 240 Hz clock interrupt request, the program determines whether the transducer reference voltage is at ground level, meaning that the switch in the external transducer is off. When the switch is first turned on, increasing the reference level to about 2.8V, the monitor program enters the main data collection, preprocessing and recording routines. In this state, each 240 Hz interrupt request causes the program to read the signal and reference value of the transducer from the ADC, compute the difference (ie, the signed amplitude of the signal relative to its resting level), and sum this value. After each four interrupts, the sum of the preceding four conversions is entered into a differencing function. The sum of every four conversions at 240 Hz represents the 60 Hz digital filter developed in Experiment 5. The differencing function is equivalent to algorithm 5 described in Experiment 5. Specifically, the input to the differencing function is multiplied by 10, and the sum of the preceding and following 5 inputs are subtracted from this product. Then the absolute values of 120 such transformations are summed to form a 2-second epoch score. The largest epoch score for each minute is stored sequentially in memory.

Also, each minute, a voltage indicating the illumination level from the photocell is digitized and stored, and a time code is signaled through the LED. The time code permits external monitoring of the system operation and coordination with polygraphic recording. The monitor program fills 448 12-bit memory locations, leaving 5696 locations available for data storage. This allows us to store two 12-bit data values (activity and illumination) each minute for 47 hours and 28 minutes. Since 4-bit resolution might be adequate, up to 6 times this duration or about 12 days data might theoretically be stored, were the illumination data sacrificed and were battery changes feasible.

SLEEP RECOGNITION SOFTWARE

Having measured, processed, and recorded activity data, the remaining task for a sleep monitor system is to interpret the activity information as sleep or wake. In our prototype system, the contents of the recorder memory were sent to an Apple II microcomputer, which stored the data on discs. The data were then transferred to a more powerful Hewlett-Packard 1000 system for automatic sleep scoring. Automatic scoring used the sleep recognition algorithm which was

developed in Experiment 6.

Fourteen digital activity records totalling 12,739 minutes were collected from healthy male college students participating in experiments in which EEG, EOG, and EMG were being recorded. Time codes from the monitor via the LED were written on the polygraph record each minute. Since the physical parameters of the transducer-amplifier-ADC differed somewhat from those used previously, the parameters of the sleep recognition algorithm were recalibrated using the method described in Experiment 6, the computer deriving the optimal sleep recognition algorithm. For these recordings, the parameters producing the best agreement with EEG scoring were:

D = .036 x [.07y((i-5) + .08y(i-4) + .10y(i-3) + .11y(i-2) + .12y(i-1) + .14y(i) + .09y(i+1) + .09y(i+2) + .09y(i+3) + .10y(i+4)]

where y(i) represents the greatest 2-second epoch score in minute i, etc, and D>1.0 is scored wake. The results of this primary scoring algorithm were then subjected to the rescoring algorithm described in Experiment 6, in which corrections for sleep onset latency and brief periods of inactivity were made. Results are summarized in Table 3. Overall minute-by-minute agreement was 93.88%. The activity estimate of % Sleep exceeded the EEG estimate by an average of 1.03%, corresponding to 9.35 minutes per record, or 14.81 minutes per 24-hour day. The correlation coefficient between the two estimates of % Sleep was r=0.9724.

An additional 14 records totalling 22,514 minutes were scored prospectively using the parameters found optimal in the calibration phase. Twelve subjects were healthy male college students, while two (subjects 12 and 13) were patients with sleep complaints undergoing clinical sleep evaluations. The patient records were included to test the system on highly disturbed sleep records. As seen in Table 4, overall agreement with EEG scoring was 93.04%. The activity estimate of % Sleep exceeded the EEG estimate by an average of 1.20%, corresponding to 19.29 minutes per record or 17.27 minutes per 24-hour day. The correlation coefficient was r=0.9692.

Figure 10 illustrates the scoring of one of the activity records in the prospective series (record 6). The 48-hour record is divided into four 12-hour segments in this figure, and hours are marked along the bottom of each segment. The top trace in each segment plots the maximal activity value for each minute. Below this is the result of EEG scoring of the record (solid indicates wake, light indicates sleep), and below this is the result of automatic scoring of the activity data. The lowest trace in each segment illustrates the agreement between EEG and activity scoring (solid indicates agreement).

ASSESSMENT OF VALIDATION

The results of our prospective validation of the automatic scoring system are comparable to those obtained when the parameters were optimized retrospectively. These results indicate the performance we would expect in an operational application. The wearable sleep monitor system measures sleep with approximately the same precision as the experimental laboratory system developed in Experiments 5 and 6, and visual scoring of analog activity records in

Table 3.

Record durations, proportion of the record for which EEG and automatically scored activity scores agree, proportion of the record scored as sleep by each technique, and correlation coefficient between the proportion scored sleep by the two methods. Activity data was collected with the wearable monitor. Scoring parameters were derived retrospectively by analysis of these data.

Record	Recording Duration (minutes)	% Agreement	% Sleep (EEG)	% Sleep (<u>Activity</u>)
1	357	91.04	85.71	99.68
2	2838	95.81	18.68	19.84
3	1463	98.63	24.27	24.13
4	2839	98.20	20.61	21.91
5	452	90.27	97.35	92.04
6	335	95.22	97.61	97.01
7	456	80.70	88.16	75.88
8	491	90.63	91.85	87.37
9	493	90.06	90.06	82.56
10	478	91.84	96.03	95.82
11	474	85.02	82.70	96.41
12	494	83.40	81.58	92.51
13	1091	93.13	35.47	36.66
14	478	89.96	89.33	98.12
	Total 12739	93.88	46.38	47.41

r=0.9724

Figure 10. A 48-hour activity record collected with the wearable monitor (divided into four 12-hour segments). Hours are marked along the bottom of each segment. The top trace in each segment plots the digital activity value for each minute. Below this is the EEG score (wake is solid, sleep is light) and below this the automatic activity score. The lower trace illustrates agreement between EEG and activity scoring (solid is agreement).




Table 4.

Record durations, proportion of the record for which EEG and automatically scored activity scores agree, proportion of the record scored as sleep by each technique, and correlation coefficient between the proportion scored sleep by the two methods. Activity data was collected with the wearable monitor. Scoring was prospective, using parameters derived from previous data.

Record	Recording Duration (minutes)	% Agreement	% Sleep (EEG)	% Sleep (Activity)
1	469	95.31	94.46	93.60
2	459	76.25	92.37	71.68
3	495	93.74	93.33	93.94
4	1302	97.00	26.57	26.80
5	2838	96.16	24.07	24.95
6	2838	97.08	25.19	26.71
7	2838	87.14	10.47	23.33
8	2593	93.21	25.92	25.38
9	2675	93.61	36.11	35.33
10	1335	93.41	34.23	29.21
11	2261	95.98	40.82	41.84
12	625	80.96	66.88	74.40
13	484	85.12	84.50	89.05
14	1302	92.93	38.56 .	33.95
• 7				
Total	22514	93.04	34.28	35.48

Experiment 1. The overall accuracy is more than adequate for practical utilization.

Throughout the research and development phase of this project, our goal has been to maximize the minute-by-minute agreement between EEG and activity scoring of sleep. We feel that the 90%+ agreement we have obtained in every phase of our research is very good, and it is unlikely that any alternative activity-based scoring system could reach substantially higher agreement. One reason for this is that EEG scoring is itself an estimate. Although rules have been formulated to assist in the standardization of EEG sleep scoring, scoring is nevertheless a matter of interpretation. Two raters scoring the same EEG record seldom agree 100% of the time (Mullaney et al, 1980). Some of the activity system's discrepancies with EEG estimates are certainly due to unreliability in the standard EEG scoring methods.

A second factor reducing the agreement of EEG and activity scoring in the present studies involves errors of resolution which occur in reducing the continuous processes of sleep and wake to discrete epochs. EEG recordings are usually scored in 20-30 second epochs. In our activity scoring system, the epoch was one minute. Our criterion for assigning a score of wake or sleep to a minute required that at least 50% (30 sec) be in that state. The EEG and activity methods had somewhat different responses to minutes containing brief awakenings. A final problem involved the inescapable fact that people sometimes do not move at all during times when they are awake.

Most investigators would probably assume that EEG scoring corresponds more accurately to behavioral and functional sleep than does activity scoring. While we have no objective comparisons of EEG and activity scoring as indicators of behavioral sleep, we suspect that where the EEG and activity scores were discrepant, on some occasions, behavioral criteria would favor the activity scoring. Nevertheless, reliability data indicates that the agreement between EEG and activity estimates of TST are only slightly less than agreement of two scorers rating the same EEG record. A more important issue is the well-known subject-to-subject and night-to-night variability which is observed in sleep recordings. In general, the increased quantity of data which could be gained for a given expense by automated wrist activity scoring would more than compensate in increased precision for any loss of accuracy as compared to EEG scoring.

Even if the EEG method is slightly more accurate in estimating sleep time, other factors such as simplicity and practicality make activity scoring more desirable for many purposes. The activity transducer is more comfortable and convenient for the subject, who can sleep outside the laboratory if desired, free from electrodes. The adaptability of the activity scoring method to naturalistic settings also increases its value as a measurement approach.

Although our current prototype activity scoring device is rather large and fragile for field conditions, we are certain that a device of equivalent capability could be miniaturized so that it could be worn entirely on the wrist-A wrist-mounted device could be built for military specifications of durability.

HARDWARE SPECIFICATIONS OF A WRIST-WEARABLE SLEEP MONITOR

Essentially all of the design, software, and procedural requirements have now been developed so that a miniaturized wrist-wearable sleep monitor can be built suitable for military settings. Although the actual construction, testing, and deployment of such a device is beyond the scope of this contract, we present below technical requirements for such a device. We believe a fully wrist-wearable device can be built with sufficient durability and reliability for military settings, and an appropriate read-out device can be constructed. Either current flat-pack or chip carrier technology would be adequate for miniaturization, the major problem being conservation of power and battery life.

The miniaturized sleep monitor should contain the following components or their equivalent:

- 1. Low-power microprocessor, processing at least 100,000 instructions per second, 8 bit word width or more.
- 2. Crystal-controlled hardware interrupts at 240 Hz or some multiple (1200 Hz recommended).
- 3. 8-K bytes of RAM memory or more.
- 4. A piezo-ceramic activity transducer as described in Figure 9, with appropriate buffer amplifier circuit and level shifter.
- 5. A single-channel A/D converter sensitive to positive and negative voltages or 2-channel positive-only A/D converter, 8-bit resolution or better and conversion rate faster than 240 Hz.
- 6. A facility for loading programs down into the microprocessor and restarting the microprocessor, either using a ROM controlled peripheral interface element or using direct control of the memory bus to load memory, eg, communication through DMA.
- 7. A facility for halting processor operations and unloading memory to an external controller, either thrugh a peripheral interface element controlled by ROM or through DMA.
- 8. Total system power consumption and battery supply must allow operation for at least 48 hours with facility for rapid battery change. Power consumption of less than 3 milliamps is recomended. Longer battery life is highly desirable.

The following design features are recommended but not required:

- 1. Program down-loading and data dump could be accomplished through optical coupling, using a self-clocking system. Two LEDs and two photocells (LED receivers) should be sufficient to implement 2-way communication. The LEDs could be controlled by I/O lines, for example, controlled through a peripheral interface element, and the software implementing this control could be housed in ROM. The photocells should have an opaque cover when not in use.
- 2. The unit should have a waterproof seal for all elements except the battery element should. The battery compartment should also have a waterproof gasket.
- 3. Parallel battery leads should be provided so that a fresh or recharged battery can be connected before a partially discharged battery is disconnected, to implement continuous operation.

4. Design of the circuit with an 8-channel multiplex analog digital converter would be useful to allow for future design expansions and implementations.

5. Inclusion of an ambient-light photocell with a window through the case, amplifier circuit, and ADC channel would be useful for sensing subject light exposure, eg, when the subject goes to bed. A sensitivity range of 1-5,000 lux or up to 50,000 lux if possible is desired.

6. If an 8-channel ADC is implemented, 3 thermistor bridges with amplifiers could be included. An external plug could allow one external thermistor to be connected. Thermistors on the top and bottom of the case could be used to detect when the instrument is taken off the wearer.

7. If the microprocessor had a low-power halt mode and could be restarted by interrupts, this capability could be used to reduce power consumption.

8. Ideally, the case should be cigarette-pack size or smaller, eg,
2x5x8 cm maximum. In our experience, a very broad leather or plastic
wrist band as wide as the case is most confortable.

9. A 4-digit LCD display controlled by the microprocessor would be extremely useful for clock and warning functions.

10. If the monitor will communicate with its external controller through DMA, the external interface must be designed. We would favor building an interface to a field-portable microprocessor such as an Apple II computer. The complete sleep scoring, readout, and display functions could be performed by a microprocessor comparable to an Apple II (16K byte 6502 processor with keyboard and display).

SOFTWARE SPECIFICATION

The logical steps necessary for sleep scoring are outlined in the preceding text and are specified exactly in the accompanying program (Appendix 1).

The software will need to be adapted for the exact configuration of wrist-wearable sleep monitor which is built, and a scaling parameter will need to be determined empirically. Operational requirements for data reporting and display will also require certain modifications of our current programs.

Now that the basic principles have been established, the hardware/software requirements for a readout device are not great. A portable microprocessor comparable to a 16K Apple II with one 5-1/4 inch disc would be sufficient for field requirements, and could easily be programmed once the report configuration desired is specified.

SUMMARY

Experiment 1 tested the idea that sleep/wake can be recognized from activity data. One hundred and two analog recordings of REG-EOG-EMG and wrist activity were collected on cassette taps. Polygraph write-outs of the EEG-EOG-EMG channels and the activity channel were scored for sleep by

independent raters. The raters agreed 94.5% of the time, indicating that wrist activity data provides nearly the same information as electrophysiological measures for recognizing sleep.

The adequacy of the piezo-ceramic activity transducer used in Experiment 1, and its orientation and placement were investigated in Experiments 2 through 4. In Experiment 2, the transducer was compared with a tilt-switch transducer and a sensitive accelerometer and was found to be superior. Experiment 3 found the transducer to be essentially omni-directional, thus making its orientation unimportant. The suitability of the wrists for placement of the transducer was confirmed in Experiment 4.

While more efficient than EEG scoring, the analog activity recordings collected in the preceding experiments required much time for playback and interpretation by a trained scorer. A practical sleep monitor system would quantify the activity data and make it available to a computer for automatic scoring. Experiment 5 derived an algorithm for reducing two seconds of activity to a single value. The optimal algorithm also effectively eliminated any 60 Hz interference.

Experiment 6 investigated methods of interpreting the digital activity value for each two-second data epoch as sleep or wake. It was found that the largest two-second epoch score in each minute best represented activity in that minute. A given minute was scored wake if a weighted sum of this score for the given minute and for the immediately preceding and following minutes exceeded a threshold. For the 17 activity records with which these rules were dervied, automatic scoring agreed with hand EEG scoring 94.46% of the 11,002 minutes. A further three records scored prospectively with the same parameters agreed 96.02% of the time with EEG scoring.

A prototype sleep monitor system was developed. The system consists of a wrist-mounted piezo-ceramic transducer, microprocessor-controlled data compression and storage device worn on the belt, and sleep recognition software in an external computer. The unit is capable of recording wrist activity digitally for several days on a single battery charge while subjects go about their normal routine. Fourteen subjects, records were used to retrospectively derive the optimal sleep recognition algorithm for this system. Of the 12,739 minutes recorded, 93.88% were scored the same by automatic scoring and by hand EEG scoring. An additional group of 14 records were scored prospectively using the same parameters, and 93.04% of the minutes were scored the same by automatic and by hand EEG techniques. The proportion of the record scored as sleep automatically from activity and from EEG differed by 1.03% and 1.20% in the two groups, and the correlation coefficients between the two estimates of sleep across subjects were .9724 and .9692.

The electronic components of this sleep monitor system could be miniaturized with current technology into a package worn entirely on the wrist. A miniature wrist-mounted sleep monitor system would make large-scale sleep measurement feasible under field conditions.

We are extremely pleased with the success of this technical development and hope our designs can be operationally applied.

ACKNOWLEDGMENT

This research was supported in part by the Veterans Administration, by ONR NO0014-79-C-0317, and by Research Scientist Development Award 1 KO1 MH00117 (to DFK) from the National Institute of Mental Health. William Mason, M.A. and Nobuyuki Okudaira, M.D., Ph.D. assisted in the research.

REFERENCES

- Carskadon, M.A., Dement, W.C., Mitler, M.M., et al. Self-reports versus sleep laboratory findings in 122 drug-free subjects with complaints of chronic insomnia. Am. J. Psychiatry, 133: 1382-88, 1976.
- Foster, F.G., Kupfer, D., and McPsrtland, R. Telemetric mobility recording in assessment of day and night cycles. Sleep Research, 1: 171, 1972a.
- Foster, F.G., Kupfer, D., Weiss, G., Lipponen, V., McPartland, R., and Delgado, J. Mobility recording and cycle research in neuropsychiatry.

 J. Interdiscip. Cycle Res., 3(1): 61-72, 1972b.
- Kripke, D.F., Messin, S., Mullaney, D.J., and Wyborney, V.G. Wrist actigraph measures sleep. In: Chase et al. (eds), <u>Sleep Research</u>. Los Angeles: Brain Information Service/Brain Research Institute, University of California, 1977, p. 210.
- Kripke, D.F., Mullaney, D.J., Messin, S., and Wyborney, V.G. Wrist actigraph measures of sleep and rhythms. <u>Electroenceph. Clin. Neurophysiol.</u>, 44: 674-676, 1978.
- Mullaney, D.J., Kripke, D.F., and Messin, S. Wrist-actigraphic estimation of sleep time. Sleep, 3(1): 83-92, 1980 (supported in part by the DAMD).

APPENDIX 1

Monitor program used to control data collection, processing, and storage for wearable sleep monitor.

```
18881 -- 28
                      /ACTICORD
5852
2223
                       /RECORDS OUTPUT OF ACTIGRAPH TRANSDUCER AND PHOTOCELL
2224
                       /RESOLUTION=1 MINUTE, DURATION=47HR,28MIN
8883
                       /TIME CODE EACH MINUTE THRU L.E.D.
3336
3337
                       /AUTHOR: JOHN WEBSTER
2222
8889
8818
5511
SS12
SS13
                      STORE - JMS I STORE
5514
8815 ·
                      /PAGE #
##15
8817
8818
                       •#
8819
8828 888 8888
                       RTNADR,
                                8
                                                 /INTRPT RETURN ADDR
                       INTRPT, JMP
##21 ##1 5#26
                                       INTSRV
                                                 /INTRPT PNTR
8822
##23
8824
                      /AUTOINDEX REGISTERS
5525
                       *13
8826 813 8888
                       A13.
                                8
                                                 /RESERVED FOR A/D. . .
8827 814 8888
                       A14.
                                8
                                                 /BUFFER POINTER
8828 815 8888
                                                 /RESERVED FOR DIFF. . .
                       A15.
                                8
                              s
##29 #16 ####
                       A16,
                                                 /SCORE POINTERS
                      DATAP. #
8838 817 B888
                                                 /DATA STORE POINTER
SS31
5532
8833
                       /INTERRUPT SERVICE ROUTINE
SE34
                       *26
##35 #26 31#6
                       INTSRV.
                                DCA
                                        ACSAVE /SAVE REGISTERS
£#36 #27 7##4
                                RAL
8837 838 3187
                                DCA
                                        LKSAVE
##38 #31
                                RCRA
           6584
8839 832 3118
                              DCA
                                        CRASAV
                                        GETACT /GET ACTIVITY SCORE
8848 833 4844
                                JMS
8941 834 1118
                                        CRASAV /RESTORE REGISTERS
                                TAD
##42 #35
                                WCRA _
           6585
8843 836 6887
                                CAF
##44 #37 11#7
                                        LKSAVE
                                TAD
8845 848 7818
                                RAR
2845 841 1186
                                TAD
                                        ACSAVE
8847 842 6881
                                ION
                                                  /ENABLE INTRPT
8848 843 5488
                                JMP I RTNADR
3549
8858
8851
                       /ACTIVITY SCORE SET-UP
8852 844
         7482
                       GETACT, HLT
8853 845 1112
                                TAD
                                       AWORD /SET UP PIE FOR A /D
##54 #46 65#$
                                WCRA
                                                  /START WITH CH.1 (BIAS)
E855 847
           7251
                                CLA IAC
8856 858
           4566
                                JMS
                                        AZD
8857 851
8858 852
           7541
                                CIA
                                        BIAS /NEGATE AND SAVE
           3121
                                DCA
2259 252 4266
                                JMS
                                        AZD
                                                /THEN CH.# (ACT)
3868 854
                                        BIAS -
          1121
                                TAD
                                                 /FIND DIFFERENCE
                                DCA I A13
TAD A13
##61 #55
                                                 /PACK INTO BUFFER
         3413
5562 556
          1513
                                                 /RESET BUFFER PHTR. . .
8863 857
                                        STKEND /IF AT END (EOB)
           1114
                                TAD
3364 363
                                SZA CLA 36
          7648
```

```
B865 B61
            5#54
                                  JMP
                                           .+3
 3366 A62
            1111
                                  TAD
                                           STEREC
 3857 863
            3513
                                  DCA
                                           A13
 #268 #64
            5444
                                  JMP I
                                          READAD
 0069
 3275
 $871
                        /ANALOG TO DIGITAL CONVERSION
 EJ72
                        *66
 3873 B56
          7452
                        AZD.
                                  HLT
 8874 B67
            65#1
                                 VRITE1
                                                    /CHAN. SELECT
 3875 878
           7388
                                 CLA CLL
 9976 B71
           1115
                                 TAD
                                          MCNT
                                                    /SET UP BACKUP CMTR
 5877 $72
           3116
                                 DCA
                                           MCNTR
 8878 873
           6511
                                 WRITE2
                                                    /START CONVERSION
 8879 874
           6513
                                 SKIP2
                                                    /DONE (EOC)?
 8888 875
           7418
                                 SKP
 ##81 #76
           5873
                                 JMP
                                           DONE
 BS82 877
           2116
                                 ISZ
                                          MCNTR
                                                    /NO, INCR. BACKUP
 9983 183
          5074
                                 JMP
                                          . - 4
SE84 181
          6500
                        DONE.
                                 READ1
                                                    /YES
$$85 152
          $117
                                 AND
                                          X377
                                                   /MASK GARBAGE &
BB86 183
          5466
                                 JMP I
                                          AZD
                                                   /EXIT
2887
8888
3389
                        /CONSTANTS AND VARIABLES
8898 184 8238
                                          /ACCUMULATES MAXIMAL EPOCH SCORE
                        MAXEP.
                               8
2891 185 8288
                        LITE,
                                8
                                         /HOLDS LITE LEVEL VALUE
                        ACSAVE, 8
                                        /SAVES AC DURING INTRPT
/SAVES LINK DURING INTRPT
3092 186 8988
9993 187 8858
                        LKSAVE. #
8894 118
                        CRASAV, #
STKBEG, 547
          2262
                                        /SAVES CRA DURING INTRPT
8895 111
          #547
                                         /POINTS TO ACT SCORE BUFFER
$896 112
          2882
                        ERROR, #
##97 113 #2#1
                        AWORD, 281
                                        /INITIAL CRA VALUE
##98 114
          7221
                        STKEND, 7221
                                        /-END OF ACT SCORE BUFFER
                       MCNT, 7766
MCNTR, #
X377, 377
M6, -6
8899 115
          7766
                                         /TIME DELAY FOR A/D
$198 116 BBBB
                                          /COUNTER FOR TIME DELAY
B181 117
          #377
                                          /USED TO MASK A/D CONV
$182 128
          7772
$183 121
          8388
                        BIAS, #
                                          /HOLDS TRANS. BIAS VALUE
S184 122
          4228
                        BWORD, 4228
                                          /INITIAL CRB VALUE
#195 123
          8677
                        DATBEG. 677
                                          /POINTS TO START OF RECORDER MEM.
S186 124
          8938
                                          /INDICATES MEMORY FIELD (# OR 1)
                        FLDFLG, #
#187 125
          5005
                       MINCHT, S
                                          /MINUTE COUNT
J188 126
          8888
                       TCODE.
                                8
                                          /TIME CODE
Ø189 127
          7775
                               -3
                       МЗ,
3118 138
          8388
                       GRPS,
                                          /COUNTER FOR TIME CODE
#111 131
          5555
                       OCNTR, 5
Ø112 132
          7773
                       M5, -5
$113 133 BEES
                       TCNTR. #
3114 134 $888
                        CODE1, #
                                          /TIME CODE '8' OR '1'
$115 135 BEES
                       TFLAG, # -3#
                                          /TIME CODE FLAG
3116 136 7742
3117 137 8858
                       MINDIV.
                                         /MINUTE DIVIDER
3113 148
          2225
                       EPSUM, #
                                          /ACCUMULATES EPOCH SUM
B119 141 2508
                       DFSUM, #
                                          /ACCUMULATES DIFSCR SUMS
5128 142
          7618
                       M128. -128
3121 143 SSES
                       EPDIV.
                                         /EPOCH DIVIDER
J122 144 BS38
                       SUM4.
                               .
                                         /ACCUMULATES SUM TO 4
#123 145
         7774
                       M4,
                               -4
$124 146 $888
                               •
                       CNT4.
                                         /COUNTS 4
                             177
J125 147 #177
                       X177.
                                         /MASK FOR DIFSCR
3126 158
                              7788
                       X7755,
          7755
S127 151
          2358
                       DATURD, #
                                          /HOLDS DATA PRIOR TO STORE
#128 152
          3688
                       CODETS,
                                CODET
                                         /POINTS TO TIME CODE RTN
```

#129	8488 8888 1888 2881 7776 7766 7688 4888 8888 8888 8888 8888	DFSC. CRATMP. FLD1. MEMLIN. M2. M18. T288. 04888. LATEST. BUFLOC. BUFN. MINUS. OFLOW. PRDCT. STORE. DIFOUT. BUMP. M11. M8. M1. X77.	DIFSCR # 1888 2881 -2 -18 7688 4888 # 557 # # # \$ \$ \$ \$ \$ \$ \$ \$ \$ \$ \$ \$ \$ \$ \$ \$	/TEMP CR. /WCRA TO /- END OF /WCRA TO /HOLDS NO /POINTS /BUFFER	GET FIELD 1 F FIELD WAIT DST RECENT SUM4 FO DIFSCR BUFFER
### ### ### ### ### ### ### ### ### ##	#### 1888 2881 7776 7766 7688 4888 8888 8888 8888 8888	CRATHP. FLD1, MEMLIH, M2, M18, T288, O4888, LATEST, BUFLOC, BUFN, MINUS, OFLOW, PRDCT, STORE, DIFOUT, BUMP, M11, M8, M1,	# 1### 2##1 -2 -1# 76## # # # # # # # # # # # # # # # # #	/TEMP CR. /WCRA TO /- END OF /WCRA TO /HOLDS NO /POINTS /BUFFER	A STORE GET FIELD 1 F FIELD WAIT DST RECENT SUM4 TO DIFSCR BUFFER LOC
### ### ### ### ### ### ### ### ### ##	1888 2881 7776 7766 7688 4888 8888 8888 8888 8888	FLD1, MEMLIM, M2, M18, T288, O4888, LATEST, BUFLOC, BUFN, MINUS, OFLOW, PRDCT, STORE, DIFOUT, BUMP, M11, M8, M1,	1888 2881 -2 -18 7688 4888 8 557 8 8 5 STORES 8 -11 -8	/WCRA TO /WCRA TO /HOLDS MG /POINTS TO /BUFFER MG	GET FIELD 1 F FIELD WAIT DST RECENT SUM4 FO DIFSCR BUFFER LOC
#132 186 #133 187 #134 168 #135 161 #136 162 #137 163 #138 164 #139 165 #148 166 #141 167 #142 178 #143 171 #144 172 #145 173 #146 174 #147 175 #148 176 #149 177 #15# #151 #152 #153 #154 #155 #156 #157 #158 #155 #156 #157 #158 #158 #159 #164 #162 #163 #164 #165 #167 #168 #167 #168 #167 #168 #167 #168 #167 #168 #169 #168 #167 #168 #169 #168 #161 #162 #168 #167 #168 #169 #168 #161 #162 #168 #163 #164 #165 #166 #167 #168 #167 #168 #169 #168 #168 #168 #168 #168 #168 #168 #168	2881 7776 7766 7688 4888 8888 8888 8888 8888	M2. M18. T288. 04888. LATEST. BUFLOC. BUFN. MINUS. OFLOW. PRDCT. STORE. DIFOUT. BUMP. M11. M8. M1.	-2 -18 7688 4888 8 557 8 8 STORES 8 -11 -8	/WCRA TO /HOLDS MG /POINTS /BUFFER G	WAIT DST RECENT SUM4 TO DIFSCR BUFFER LOC
#133 157 #134 16# #135 161 #136 162 #137 163 #138 164 #139 165 #148 166 #141 167 #142 17# #144 172 #145 173 #146 174 #147 175 #148 176 #149 177 #15# #151 #152 #153 #154 #155 #156 #157 #158 #158 #159 #164 #155 #164 #165 #167 #168 #167 #168 #167 #168 #167 #168 #167 #168 #167 #168 #167 #168 #169 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #163 #164 #165 #166 #167 #168 #167 #168 #169 #168 #168 #169 #168 #168 #168 #168 #168 #168 #168 #168	7776 7786 7688 4888 8888 8557 8888 8888 8888 8648 8888 7765 7778	M18. T288. 04888. LATEST. BUFLOC. BUFN. MINUS. OFLOW. PRDCT. STORE. DIFOUT. BUMP. M11. M8. M1.	-18 7688 4888 8 557 8 8 8 STORES 8 -11 -8	/HOLDS MI /POINTS T /BUFFER (OST RECENT SUM4 FO DIFSCR BUFFER LOC
#135 161 #136 162 #137 163 #138 164 #139 165 #148 166 #141 167 #142 17# #143 171 #144 172 #145 173 #146 174 #147 175 #148 176 #149 177 #15# #151 #152 #153 #154 #155 #156 #157 #158 #158 #159 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #167 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #163 #164 #165 #166 #167 #168 #169 #168 #169 #168 #169 #169 #168 #169 #178 #178 #178 #178 #178 #178 #178 #178	7786 7688 4888 8888 8557 8888 8888 8888 8648 8888 8888 7765 7778	T288. 04888. LATEST. BUFLOC. BUFN. MINUS. OFLOW. PRDCT. STORE. DIFOUT. BUMP. M11. M8. M1.	76## #### # 557 # # # STORES # -11 -8 -1	/HOLDS MI /POINTS T /BUFFER (OST RECENT SUM4 FO DIFSCR BUFFER LOC
#136 162 #137 163 #138 164 #139 165 #148 166 #141 167 #142 17# #143 171 #144 172 #145 173 #146 174 #147 175 #148 176 #149 177 #15# #151 #152 #153 #154 #155 #156 #157 #158 #158 #159 #168 #161 #162 #168 #163 #164 #165 #166 #167 #168 #168 #169 #168 #169 #169 #169 #169 #169 #169 #169 #169	4888 8888 8557 8888 8888 8888 8648 8888 8888 7765 7778	O4###. LATEST. BUFLOC. BUFN. MINUS. OFLOW. PRDCT. STORE. DIFOUT. BUMP. M11. M8. M1.	######################################	/HOLDS MI /POINTS T /BUFFER (OST RECENT SUM4 FO DIFSCR BUFFER LOC
#137 163 #138 164 #138 165 #148 166 #141 167 #142 17# #143 171 #144 172 #145 173 #146 174 #147 175 #148 176 #149 177 #15# #153 #154 #155 #156 #157 #158 #158 #159 #168 #161 #162 #168 #163 #164 #165 #166 #167 #168 #168 #169 #168 #169 #169 #169 #169 #169 #169 #169 #169	### #557 ### ### ### ### #648 ### ### ### ### ### ### ### ### ### #	LATEST, BUFLOC, BUFN, MINUS, OFLOW, PRDCT, STORE, DIFOUT, BUMP, M11, M8, M1,	# 557 # # # # # # # # # # # # # # # # #	/HOLDS MI /POINTS T /BUFFER (OST RECENT SUM4 FO DIFSCR BUFFER LOC
#13# 164 #139 165 #14# 166 #141 167 #142 17# #143 171 #144 172 #145 173 #146 174 #147 175 #148 176 #149 177 #15# #151 #152 #153 #154 #155 #156 #157 #158 #158 #159 #164 #162 #163 #164 #165 #167 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #163 #164 #165 #166 #167 #168 #169 #168 #169 #169 #169 #169 #169 #169 #169 #169	#557 ### ### ### #8## #64# ### ### ### ### ### ### ### ###	BUFLOC, BUFN, MINUS, OFLOW, PRDCT, STORE, DIFOUT, BUMP, M11, M8, M1,	557 8 8 8 STORES 8 -11 -8 -1	/POINTS 1 /BUFFER 1	TO DIFSCR BUFFER LOC
### 165 ### 166 ### 166 ### 166 ### 166 ### 167 ### 167 ### 167 ### 173 ### 173 ### 173 ### 174 ### 175 ### 175 ### 176 ### 177 ### 175 ### 176 ### 177 ### 175 #### 176 #### 177 #### 175 #### 176 #### 177 #### 175 #### 176 #### 177 ##### 177 ############# ##########	8888 8888 8888 8648 8888 8888 7765 7778	BUFN, MINUS, OFLOW, PRDCT, STORE, DIFOUT, BUMP, M11, M8,	# # # # # # # # # # # # # # # # # # #	/BUFFER (LOC
#14# 166 #141 167 #142 17# #143 171 #144 172 #145 173 #146 174 #147 175 #148 176 #149 177 #15# #151 #152 #153 #154 #155 #156 #157 #158 #159 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #163 #164 #165 #166 #167 #168 #169 #168 #169 #168 #171 #168 #172 #178 #178 #178 #178 #178 #178 #178 #179 #178 #178 #178 #178 #178 #178 #178 #178	8888 8888 8648 8888 8888 7765 7778	MINUS, OFLOW, PRDCT, STORE, DIFOUT, BUMP, M11, M8, M1,	# # # # # # # # # # # # # # # # # # #		
#141 167 #142 17# #143 171 #144 172 #145 173 #146 174 #147 175 #148 176 #149 177 #15# #152 #153 #154 #155 #156 #157 #158 #158 #159 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #163 #164 #165 #164 #165 #165 #166 #167 #168 #169 #168 #169 #168 #171 #168 #178 #178 #178 #178 #178 #178 #178 #17	8888 8888 8648 8888 8888 7765 7778 7777	OFLOW, PRDCT, STORE, DIFOUT, BUMP, M11, M8, M1,	# STORES # -11 -8 -1	/POINTS	TO STORE RTH
#142 17# #143 171 #144 172 #145 173 #146 174 #147 175 #148 176 #149 177 #15# #151 #152 #153 #154 #155 #156 #157 #158 #158 #159 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #161 #162 #168 #165 #164 #165 #165 #166 #167 #168 #167 #168 #169 #169 #169 #178 #178 #178 #178 #178 #178 #178 #178	8888 8648 8888 8888 7765 7778 7777	PRDCT, STORE, DIFOUT, BUMP, M11, M8, M1,	# STORES # # -11 -8 -1	/POINTS	TO STORE RTN
### ### ### ### ### ### ### ### ### ##	8648 8888 8888 7765 7778 7777	STORE, DIFOUT, BUMP, M11, M8, M1,	STORES # -11 -8 -1	/POINTS	TO STORE RTN
### 172 ### 173 ### 174 ### 175 ### 176 ### 176 ### 176 ### 177 ### 176 ### 177 ### 177 #### 177 #### 177 #### 177 #### 177 #### 177 #### 177 ##### 177 ##########	8888 8888 7765 7778 7777	DIFOUT, BUMP, M11, M8, M1,	# -11 -8 -1	/POINTS	TO STORE RTH
### 173 ### 173 ### 174 ### 175 ### 176 ### 176 ### 177 ### 178 ### 176 ### 177 #### 177 #### 177 #### 177 #### 177 #### 177 #### 177 #### 177 #### 177 ##### 177 ##### 177 ##########	8888 7765 7778 7777	BUMP, M11, M8, M1,	# -11 -8 -1		
### ### ### ### ### ### ### ### ### ##	7765 777 8 7777	M11. M8. M1.	-11 -8 -1		
### 175 ### 176 ### 176 ### 177 ### 178 ### 177 #### 177 #### 177 #### 177 #### 177 #### 177 ##### 177 ##### 177 ##########	777 8 7777	M8. M1.	-8 -1		
### 176 ### 177 ###############################	7777	M1,	-1		
### 177 ### 177 #### 158 ####################################					
### ### ### ### ### ### ### ### ### ##	5677	x77.	77		
#151 #152 #153 #154 #155 #156 #157 #158 #159 #168 #161 #162 #163 #164 #165 #166 2## #165 #166 2## #172 2#6 #173 2#7 #174 21# #177 213 #178 214 #179 215 #178 214 #179 215 #178 214 #179 215 #182 22#					
#152 #153 #154 #155 #156 #157 #158 #159 #168 #161 #162 #163 #164 #165 #166 2## #165 #166 2## #172 2#6 #173 2#7 #174 21# #177 213 #178 214 #179 215 #178 214 #179 215 #178 214 #179 215 #182 22#					
#153 #154 #155 #156 #157 #158 #159 #168 #161 #162 #163 #164 #165 #166 2## #165 2#1 #168 2#2 #169 2#3 #177 2#4 #171 2#5 #172 2#6 #173 2#7 #174 21# #175 211 #176 212 #177 213 #178 214 #179 215 #188 216 #181 217 #182 22#					
#154 #155 #155 #156 #157 #158 #159 #168 #161 #162 #163 #164 #165 #166 288 #167 281 #168 282 #169 283 #178 284 #171 285 #172 286 #172 286 #173 287 #174 218 #175 211 #176 212 #177 213 #178 214 #179 215 #181 217 #182 228					
#155 #156 #157 #158 #159 #168 #161 #162 #163 #164 #165 #166 288 #167 281 #168 282 #169 283 #178 284 #171 285 #172 286 #173 287 #174 218 #175 211 #176 212 #177 213 #178 214 #179 215 #181 217 #182 228					
#156 #157 #158 #159 #168 #161 #162 #163 #164 #165 #166 288 #167 281 #168 282 #169 283 #178 284 #171 285 #172 286 #173 287 #174 218 #175 211 #176 212 #177 213 #178 214 #179 215 #188 216 #181 217 #182 228					
#157 #158 #159 #168 #161 #162 #164 #165 #166 288 #167 281 #168 282 #169 283 #178 284 #171 285 #172 286 #173 287 #174 21# #175 211 #176 212 #177 213 #178 214 #179 215 #178 216 #181 217 #182 228					
#158 #159 #168 #161 #162 #163 #164 #165 #166 288 #167 281 #168 282 #169 283 #178 284 #171 285 #172 286 #173 287 #174 218 #175 211 #176 212 #177 213 #178 214 #179 215 #181 217 #182 228					
#159 #168 #161 #162 #163 #164 #165 #166 2## #167 2#1 #168 2#2 #169 2#3 #178 2#4 #171 2#5 #172 2#6 #173 2#7 #174 21# #175 211 #176 212 #177 213 #178 214 #179 215 #179 215 #178 216 #181 217 #182 22#					
#16# #161 #162 #163 #164 #165 #166 2## #167 2#1 #168 2#2 #169 2#3 #178 2#4 #171 2#5 #172 2#6 #173 2#7 #174 21# #175 211 #176 212 #177 213 #178 214 #179 215 #181 217 #182 22#					
#161 #162 #163 #164 #165 #166 288 #167 281 #168 282 #169 283 #178 284 #171 285 #172 286 #173 287 #174 218 #175 211 #176 212 #177 213 #178 214 #179 215 #188 216 #181 217					
#162 #163 #164 #165 #166 288 #167 281 #168 282 #169 283 #178 284 #171 285 #172 286 #173 287 #174 218 #175 211 #176 212 #177 213 #178 214 #179 215 #188 216 #181 217		/PAGE 1			
\$153	•	7.4.55			
#164 #165 #166 288 #167 281 #168 282 #169 283 #178 284 #171 285 #172 286 #173 287 #174 218 #175 211 #176 212 #177 213 #178 214 #179 215 #188 216 #181 217		•			
### ### ### ### ### ### ### ### ### ##		/START-U	P		
### ### ### ### ### ### ### ### ### ##		*255			
### ### ### ### ### ### ### ### ### ##	6887	START,	CAF		
#169 283 #178 284 #171 285 #172 286 #173 287 #174 21# #175 211 #176 212 #177 213 #178 214 #179 215 #188 216 #181 217 #182 22#	1113		TAD	AWORD	SET UP PIE
### 284 ### 285 ### 286 ### 286 ### 286 ### 287 ### 218 ### 211 ### 212 ### 213 ### 214 ### 215 ### 216 ### 217 ### 218 ### 216 ### 218 ### 217 ### 218 ### 218 ### 218 ### 218 ### 218 ### 218	65#5		WCRA		
### ### ### ### ### ### ### ### ### ##	7388		CLA CLL		
#172 286 #173 287 #174 218 #175 211 #176 212 #177 213 #178 214 #179 215 #188 216 #181 217 #182 228	1122		TAD	BWORD	
#173 2#7 #174 21# #175 211 #176 212 #177 213 #178 214 #179 215 #188 216 #181 217 #182 22#	6515	•	WCRB		
\$174 218 \$175 211 \$176 212 \$177 213 \$178 214 \$179 215 \$188 216 \$181 217 \$182 228	7355		CLA CLL		
\$175 211 \$176 212 \$177 213 \$178 214 \$179 215 \$188 216 \$181 217 \$182 228	1123		TAD	DATBEG	/ZERO MEMORY
\$176 212 \$177 213 \$178 214 \$179 215 \$188 216 \$181 217 \$182 228	3817		DCA	DATAP	
\$177 213 \$178 214 \$179 215 \$188 216 \$181 217 \$182 228	3151		DCA	DATWRD	
#178 214 #3179 215 #318# 216 #3181 217 #3182 22#	4571		STORE		•
3179 215 3188 216 8181 217 8182 228	7418		SKP	_	
5185 216 5181 217 5182 225			JMP	2	MILTO PAR 1441
\$181 217 \$182 22\$	5212	WAIT,	CLA IAC	400	/WAIT FOR 'ON'
\$182 22#	7251		JMS	AZD	\$.
	72 51 4 868		TAD	M1 <i>88</i>	ATE CUITCH ONS
miss ZZl	7281 4868 1158		SMA CLA	60	/IS SWITCH ON?
8184 222	7281 4868 1158 7788		JMP RCRA	60	/NO, WAIT
#185 223	7281 4868 1158 7788 5226			04888	INV. WALL
#186 224	7281 4868 1158 7788 5226 6584		TAD	V-222	
\$187 225	7281 4868 1158 7788 5226 6584 1162		MCDA	WAIT	
#188 226	7281 4868 1158 7788 5226 6584 1162 6585		WCRA JMP	WC-1	/YES, INIT 4 GO
#189 227	7281 4858 1158 7788 5226 6584 1162 6585 5215		JMP	AUNER -	* 1541 1011 # 44
Ø198 238	7281 4868 1158 7788 5226 6584 1162 6585 5218	60,	JMP TAD	AWORD -	
\$191 231	7281 4868 1158 7788 5226 6584 1162 6585 5215 1113	60 ,	JMP TAD WCRA	AWORD -	
#192 232	7281 4868 1158 7788 5226 6584 1162 6585 5218 1113 6585 7888	60,	JMP TAD WCRA NOP	AWORD -	
	7281 4868 1158 7788 5226 6584 1162 6585 5215 1113	60,	JMP TAD WCRA	AWORD -	

#193	233	3517		DCA	DATAP	•
.4194	- 234	3124		DCA	FLDFLG	
<i>#</i> 195		3125		DCA	MINCHT	
#196		1111		TAD	. STKBEG	
#197		3513		DCA	A13	
#198	245	1111		TAD	STKBEG	
#199	241	3514		DCA	A14	
3288		1125	MINIT.	TAD	MINCHT	JENTES PAGE MANUES
5 251	243	3126	***************************************	DCA	TCODE	/ENTER EACH MINUTE
5282		1145		TAD	M4	/LOAD TIME CODE
\$283		313#		DCA	GRPS	
8284		1145		TAD		
3285		3131		DCA	M4	
#2#6		1176		TAD	OCNTR	
\$287		3133		DCA	M1	
\$2\$ 3		1176		TAD	TCNTR	
#239		3134			M1	:
#215		2135		DCA	CODE1	
#211		2125		ISZ	TFLAG	SET TIME CODE FLAG
\$212		1136		ISZ	MINCHT	
Ø213		3137		TAD	M3 <i>B</i>	
B 214	268	3184		DCA	MINDIV	
#215		3185		DCA	MAXEP	
#216	262	3145	FRACU	DCA	MINTOT	
#217	263		EPOCH,	DCA	EPSUM	/ENTER EACH EPOCH
#21B	264	3141		DCA	DFSUM	
#21B		1142		TAD	M128	
8228	265	3143		DCA	EPDIV	•
	256	3144	FOUR,	DCA	SUM4	/ENTER EACH FOUR A/D'S
3221	267	1145		TAD	M4	
#222	278	3146		DCA	CNT4	
#223	271	1135		TAD	TFLAG	
#224	272	764#		SZA CLA		
#225	273	4552		JMS 1	CODETS	
#226	274	6881		ION		
#227	275	7388	EACH,	CLA CLL		/ENTER EACH A/D
#22#	276	1813		TAD	A13	
#229	277	7841		CIA		
\$23\$	388	1514		TAD	A14	
#231	381	7448		SZA		/IS A/D BUFFER EMPTY?
<i>3</i> 23 <i>2</i>	382	5315		JMP	NEW	
Ø233	383	65#4		RCRA		/YES, ENTER WAIT STATE
8234	384	7421		MQL		
J235	3#5	7581		MQA		
2236	386	1162		TAD	04888	
<i>3</i> 237	387	65#5		WCRA	·	
2238	315	77.61		ACL		
£239	311	65#5		WCRA		
324 8	312	7888		NOP		
3241	313	7888		NOP		
#24Z	314	5275		JMP	EACH	
7243	315	7388		CLA CLL		/NO, READ NEXT VALUE
9244	316	1414		TAD I	A14	
2245	317	1144		TAD	SUM4	/AND SUM TO 4
3246	328	3144		DCA	SUM4	
3247 3249	321	1814		TAD	A14	RESET BUFFER IF AT END
£248	322	1114		TAD	STKEND	
3249	323	7648	•	SZA CLA		
D25#	324	5327	•	JMP	.+3	
7251	325	1111	·	TAD	STKREC	
Ø252	326	3814		DCA	A14	
2253	327	2146		ISZ	CNT4 -	/SUMMED 4 YET?
J254	338	5275		JMP	EACH	/NO
<i>82</i> 35	331	1144		TAD	SUM4	/YES, FIND DIFF. SCR.
3256	332	4553		JMS I	DFSCR	

	•				
. 4257	333	7188	CLL		
#258		1141	TAD	DFSUM	/ADD TO DIFSUM
#259	335	7438	SZL		
#26#	336	7248	STA	•	
\$26 1	337	3141	DCA	DFSUM	TRUNCATE IF TOO LARGE
#262	34#	2143	ISZ	EPDIV	/IS I? AN EPOCH YET?
# 263	341	5266	JMP	FOUR	/NO
3 264	342	1141 •	TAD	DFSUM	/YES, FIND MAXIMAL
#265	343	711#	CLL RAR		/EPOCH EACH MINUTE
#266	344	7841	CIA		
#267	345	1184	TAD	MAXEP	
#268	346	7788	SMA CLA		/IS THIS IT?
#269	347	5353	JMP	.+4	/NO
3273	350	1141	TAD	DF SUM	/YES, REPLACE OLD
\$271	351	7118	CLL RAR		•
\$ 272	352	3184	DCA	MAXEP	. ,
\$273	353	2137	ISZ	MINDIV	/IS IT A MINUTE YET?
<i>8</i> 274	354	5262	JMP	EPOCH	/NO
#275	355	6882	IOF		/YES. READ LITE LEVEL
	356	1363	TAD	.+5	/FROM A/D, CH.2
\$ 277		4866	JMS	A2D	
#278		3151	DCA	DATWRD	
8 279	361	6881	ION		
\$28 5	362	5367	JMP	.+5	
\$281	363	8882	2		
#282	364	7828	NOP		
#283	365	7888	NOP		
#284	366	7888	NOP		/AND 67000 IN 0500000
#285	367	4571	STORE		AND STORE IN RECORDER
\$ 285	378	7482	HLT	MAVER	/MEMORY
#287	371	1184	TAD	MAXEP Datwrd	/NOW GET MAXIMAL EPOCH /AND STORE IT IN NEXT
#288 #289	372 3 73	3151 4571	DCA Store	DAIWKD	/MEMORY LOCATION
#29#	373 374	7482	HLT		THEHORY LOCATION
#291	375	5242	JMP	MINIT	
#292	3.3	92-72	U (1)	******	
# 293					
# 294					
#295					
\$296					
#297					
#298					
8 299					
#3##					•
#3#1					
\$352					
#3#3					
3384					
#3#5					
#3#6					
#3#7					
#3#8					
2329					•
#31#					•
5311					
#312					
#313					
#314					
8315					
<i>8</i> 316				•	

in the second of the contraction of the contraction

#319 #32#

```
8321 X X
#322
                           /PAGE 2
J323
$324
J325
                           /COMPUTE DIFFERENCE SCORE
5326
#327
                           *455
                           DIFSCR.
#328 4## 74#2
                                     HLT
#329 4#1
                                     DCA
                                               LATEST
            3153
                                                       /STORE MOST RECENT SUM
#33# 4#2
                                     TAD
           1164
                                               BUFLOC
                                                         /SET UP BUFF PNTRS
#331 4#3
#332 4#4
           3516
                                     DCA
                                               A16
            1174
                                     TAD
                                               M11
8333 485
            3165
                                     DCA
                                               BUFN
#334 4#6
           3345
                                     DCA
                                               MINLO
                                                         /CLEAR DP ADDER
#335 4#7
                                     DCA
                                               MINHI
            3341
                                     CLA CLL
#336 41#
                           SUM.
                                                         /SUM; LAST 11 SUM4'S
           7355
#337 411
                                     TAD I
            1416
                                               A16
#338 412
#339 413
           7518
                                     SPA
                                     JMP
                                               NEG
           5223
8348 414
                                    TAD
                                               MINLO
           1345
                                    DCA
#341 415
           3345
                                               MINLO
#342 416
           7438
                                    SZL
#343 417
           7881
                                    IAC
8344 428
           1341
                                     TAD
                                               MINHI
#345 421
#346 422
            3341
                                     DCA
                                               MINHI
           5231
                                     JMP
                                               NEXT
8347 423
           1345
                           NEG,
                                     TAD
                                               MINLO
#348 424
                                     DCA
           3348
                                               MINLO
#349 425
           7428
                                     SNL
#35# 426
                                     STA
           7248
#351 427
                                               MINHI
           1341
                                     TAD
$352 435
$353 431
$354 432
                                     DCA
           3341
                                               MINHI
            2165
                           NEXT.
                                    ISZ
                                               BUFN
                                               SUM
           5215
                                     JMP
                                     TAD
                                               M11
#355 433
           1174
#356 434
           3165
                                     DCA
                                               BUFN
#357 435
           3342
                                     DCA
                                               PRDLO
#358 436
           3343
                                     DCA
                                               PRDHI
#359 437
           7888
                                     NOP
$36$ 44$
$361 441
           7388
                                                         /MUY CURR VAL
                           LXIS.
                                     CLA CLL
           1365
                                     TAD
                                               BUFR+5
                                                         /BY 11
#362 442
           7515
                                     SPA
#363 443
           5253
                                     JMP
                                               LNEG
#364 444
           1342
                                     TAD
                                               PRDLO
#365 445
           3342
                                    DCA
                                               PRDLO
#366 446
#367 447
#368 45#
#369 451
           7435
                                    SZL
            7851
                                    IAC
            1343
                                     TAD
                                               PRDHI
           3343
                                    DCA
                                               PRDHI
#37# 452
                                     JMP
                                               LNXT
           5261
#371 453
           1342
                                     TAD
                                               PRDLO
                           LNEG.
#372 454
           3342
                                     DCA
                                               PRDLO
8373 455
           7425
                                     SNL
#374 456
#375 457
#376 46#
           7245
                                     STA
                                     TAD
                                               PRDHI
            1343
                                     DCA
                                               PRDHI
           3343
#377 461
                                               BUFN
           2165
                           LNXT.
                                    ISZ
#378 462
           3245
                                     JMP
                                               LXIS
#379 463
           7188
                                     CLL
#38# 464
#381 465
#382 465
           1345
                                     TAD
                                               HINLO
            7841
                                     CIA
                                    DCA
                                               MINLO
            3345
8383 467
            7438
                                    SZL
8384 478
                                    IAC
            7881
```

41

Ē

	* "	••							
		471	1341		TAD		MINHI		
	5386		7841		CIA				
	#387		3341		DCA		MINHI		
	1388	474	7188		CLL				
	#389	475	1342		TAD		PRDLO		
ŀ	#39#	475	1345		TAD		MINLO	•	•
l	# 391	477	3342		DCA		PRDLO	•	
ĺ	# 392	588	7435		SZL				
	#393	551	7881		IAC				•
İ	#394	582	1343		TAD		PRDHI		
	#395	5 #3	1341		TAD		MINHI		
	#396	584	3343		DCA		PRDHI		
	#397	5#5	1164		TAD		BUFLOC		
	# 398	5.06	3816		TAD		A16	:	
	₽399	557	1164		TAD		BUFLOC	:	
	8488	515	3815		DCA		A15	•	
	\$481	511	1174		TAD		M11		
	8482	512	3165		TAD		BUFN		
	5453	513 .		SHUFL,	TAD	I	A15		
	3454	514	3173		DCA		BUMP		
	8485	515	1163		TAD		· LATEST		
	3486	516	3416		DCA	1	A16		
	8487	517	1173		TAD		BUMP		
	#458	528	3163		DCA		LATEST		
	#499	521	2165		ISZ		BUFN		
	8418	522	5313		JMP		SHUFL		
	8411	523	1342		TAD		PRDLO		
	8412	524	751#		SPA				
	8413	525	7841			•			
	8414	526	7852		BSV	-			
	8415	527	7884		RAL				
	8 416	53 #	#147		AND	_	X177		
	#417 #418	531	5688		JMP	I	DFSCR		
	8419								
	#42#			•					
	8421			/ *54#					
	#422	548	8558	MINLO,	-				
	8423	541	8888	MINHI.	•				
	8 424	542	8888	PRDLO,	•				
	#425	543	8888	PRDHI.	•				
	#426	544	8888	SCNT.	i				
	#427	545	8888	•••••	i				
	#428	546			i				
	5 429	547	8988						
	#43#	558	8888	STACK,	8				
	8431	551	8888						
	8432	552	8888						
	8 433	553	8888						
	8434	554	828 8		•				
	#435	555	8888						
	8436	356	5555						•
	#437	557	2222		8				•
	5 438	568	8888	BUFR.					
	5439	561	8888						
	8448	562	8888						
	8441	563	8888						
	5442	554	8255		•				
	8443	565	8888		•		• •		
	8444	556	8888						
	8445	567	3885		•		•		
	8446	57#	2222	•					
	8447	571	8888				42		
	BAAR	672	2222		•				

THE RESERVE AND ASSESSED ASSES

/PAGE 3

B488 /GENERATE TIME CODE (OCTAL MINUTES) 5489 *688 3495 635 7452 CODET, /ENTER EACH 4 A/D'S IF TFLAG SET HLT **3491** 551 2133 ISZ TCNTR /TIME FOR ANOTHER DIGIT? 8492 6#2 5225 JMP TOUT /NO 5493 683 116# TAD MIS 3494 654 3133 DCA TCNTR 8495 685 2131 ISZ OCNTR /YES, START OCTAL GROUP? **B496** 656 5215 JMP FLASH /NO. OUTPUT DIGIT 2497 687 1145 M4 /YES. SKIP DIGIT TAD 3498 618 3151 DCA OCNTR /4 OCTAL GROUPS YET? 5499 611 2135 GRPS ISZ 8588 612 7415 SKP 3581 613 3135 DCA TFLAG /YES, CLEAR TFLAG 8582 514 5225 JMP TOUT 3583 615 6586 FLASH, STLITE /LITE ON 3584 616 1126 TCODE TAD 8585 517 7184 CLL RAL /GET NEXT MSB 8586 628 3126 DCA TCODE J587 621 7435 SZL /BINARY 1 OR #7 jšeb 622 1127 TAD M4 /1:HOLD LITE & UNITS 8589 623 1176 TAD MI /#:HOLD LITE 1 UNIT D518 624 3134 DCA CODE 1 3511 625 1134 TOUT. TAD CODE 1 /HOLD LITE ONT J512 626 7588 SMA

: BE12	£ 627	2227						
	635	65 <i>8</i> 7 7 <i>88</i> 1			CLL			/NO. CLEAR IT
	631				JAC			/YES
	632	3134			DCA		CODE1	
	_	5688			JMP	I	. CODET	
#517								
#518								
#519				/STORE DA	ATA			
#52#				~64#			•	
#521		7452	•	STORES,	HLT			
#522		6584			RCR/	١.		/SAVE CRA
#523		3154			DCA		CRATHP	
8524		1124			TAD		FLDFLG	/WHICH FIELD?
	644	7658			SNA	CLA		
#526		5251			JMP		.+4	
#527		1155			TAD		FLD1	/IF FLD 1 CHANGE FLDS
<i>9</i> 528	647	65 <i>8</i> 5			WCRA	1		:
\$ 529		7388			CLA	CLL		•
2538		1151			TAD		DATWRD	/GET DATA
3531	652	3417			DCA	I	DATAP	
#532	653	1154			TAD		CRATMP	
#533	654	65Ø5			WCRA	•		
#534		7388			CLA	CLL		
J535		1817			TAD		DATAP	/IS MEMORY FULL?
		1156			TAD		MEML IM	
#537		765#			SNA	CLA		
Ø538		5264			JMP		.+3	
#539		2248			ISZ		STORES	
3545		5540			JMP	I	STORES	/NO, RETURN +1
2541		1124			TAD		FLDFLG	/YES, WHICH FLD?
Ø542		7448	•		SZA			
#543		5648	•		JMP	1	STORES	/FLD 1, RTRN & END
3544		2124			ISZ		FLDFLG	
£545		7248			STA			
B546		3817			DCA		DATAP	/RESET DATA PHTR TO TOP
8 547		2248			1,5Z		STORES	
	673	5648			JMP	1	STORES	/RTRN +1
#549								-
#55#								
\$ 551			4	/RECORDER	MEMOR	Y		
<i>\$</i> 552								
3553				/MAXIMAL E	POCH	ACTI	VITY SCOR	EAND
P554			4	LITE LEVE	L FOR	EAC	H MINUTE	ARE
3555				STORED IN	1 ALTE	RNAT	E MEMORY I	.ocs
Ø556				FROM LOC.	788	TO 5	777 (FLD 4	7)
Ø557 Ø558			,	AND 15555	7 TO 1	5777	(FLD 1)	•
				. =				
#\$59			4	788				
9568 9561								
3561 3562								
9562 9563			_					
~355			8	33				

END

FILMED

3-83

DTIC